OSC 97AU117 Grain Hardness Tester (Sclerometer)

- < Features >
- Measure hardness of granular good such as rice, wheat and so on.
- Detection material is stainless steel

< Specifications >

Model	OSC 97AU 117-A	OSC 97AU 117-B	OSC 97AU 117-C		
Pressure range	20kg	20kg	30kg		
Display	Dialing system	LED	LED		
Pressure method	Manual	Manual	Automatic		
max. sample size	dia.18mm	dia.18mm	dia. 20mm		
min. sample size	dia. 4mm	dia. 4mm	dia. 4mm		
Accuracy	±2%	±0.5kg	±0.5kg		
Dimension(mm)	140x100x250	330X180X100	570X330X330		
Weight	10kg	6kg	15kg		
Other function		automatic reset, automatic measuring			
		repeat, self-diagnosis			



OSC 97AU 117-A

OSC 97AU118 Fruit Hardness Tester (Sclerometer)

- < Features >
- Used for testing hardness of apple, pear, strawberry, grapefruit.
- Light weight, small volume, straight reads the type, carries conveniently,



< Specifications >

Model	OSC 97AU 118-A	OSC 97AU 118-B	OSC 97AU 118-C
Applicable fruits	Apple, pear, watermelon	Banana, fig, orange etc.	General type
Measuring range	2 ~ 15kg/cm ²	0.2 ~ 4kg/cm2	0.5 ~ 15kg/cm2
Min. scale	0.1N/cm2	0.02N/cm2	0.1N/cm2
Accuracy	±0.5N/cm2	±0.5N/cm2	±0.3 ~ ±0.5N/cm2
Measuring depth	approx.10mm	approx.10mm	approx.10mm
Tip dia.	dia. 3.5mm	dia. 3.5mm	dia.8 & dia.11mm
Dimension(mm)	approx. 160x115x50		
Weight	0.3kg		